



Adjustable Radiolucent Base Unit for Neurosurgical Headrest System

1101.021



Adjustable Radiolucent Base Unit for Neurosurgical Headrest System

- Radiolucent composition
- Direct fixation on EU rails (25x10 mm) or US rails (30 x 8mm)
- Requires one locking transitional member (1101031)
- Requires one transitional member (1101026)
- Requires swivel adaptor (3033-00 or 3033-10) to connect to head support (3034-00, 3035-00 or 3036-00)
- Ordered Upon Request - Non Stock Item

COMPATIBLE WITH

HIMAX - CMAXT - CMAXS - CMAX XRAY - GMAX - 5085 HARMONY - CMAX3 - CMAX3-DRIVE - KERNA

